


Leica Easy Dip manual slide staining system



Leica Biosystems Nussloch GmbH offers a manual, **easy-to-use low-cost** staining system. The Leica Easy Dip Slide Staining System made by Simport Plastics, Canada, is designed for customers' **manual staining requirements associated with cryosectioning and whenever the customer is in need of staining just a couple of slides at a time.**

The EasyDip™ Slide Staining System has two components: a **square staining jar and a 12-position vertical slide rack**. Jars can be loosely joined to each other laterally, therefore making sure they are kept in the same order when moved around on the lab counter. As an extra benefit, they are available in **5 different colors** to help better identifying staining fluids or reagents per jar.

The staining jar being made of resistant acetal plastic will not break like most glass jars do. It will resist to attacks by most staining agents including alcohol and xylene but not phenol.

The **EasyDip™ Slide Staining Rack will hold up to 12 microscope slides** with dimensions such as 25 x 75mm, 3 x 1" and even 26 x 76mm and with a thickness of 1.0 and 1.2mm. The slides fit into individual slots for free passage and rapid drainage of staining fluids. Since they are placed vertically in the rack and not horizontally, their writing area will not be stained by the fluid, allowing their removal without the use of forceps.



This staining rack is made of material specially formulated for **rapidly drying slides in a microwave oven**, or at temperatures up to 120°C. The lid completely covers the EasyDip™ Slide Staining Jar to minimize spill and evaporation. A handle is permanently attached to the rack for easy insertion and removal of slides without your fingers touching the solution. The base of the rack is placed in the vertical position to secure the slides in place, and is rotated sideways for allowing their easy removal.