



Silica/Ceramics particles/powder
prepared with the Leica EM TXP



Application Note

Leica EM TXP

Silica/Ceramics particles/powder

Market: Silica/Ceramic Powder manufacturer

Companies (e.g.): Gyeonqqi

Living up to Life

Leica
MICROSYSTEMS

Leica EM TXP Application No. 4/1

Silica/Ceramics particles/powder

Market: Silica/Ceramic Powder manufacturer

Companies (e.g.): Gyeonqqi, www.globalsources.com

Robert Ranner, Leica Mikrosysteme GmbH, Vienna, Austria
(Robert.Ranner@leica-microsystems.com)

Complete process time:	30 min
-------------------------------	--------

Process description:	(benchmark values for this particular sample)
-----------------------------	---

Sample embedded: No
Holder: glued on specimen stub for SEM (AFM holder insert)

Tool inserts:	r.p.m	F value	Pump	E-W speed
Diamond foils LA type: 3 μ , 1 μ , 0.5 μ for final finish Foil carrier: steel	2000	5-8	8ml	0.35mm/s



(Sample on specimen stub)

