



Leica RM2245

The Manual Microtome with High-Precision
Motorized Specimen Feed

Leica
MICROSYSTEMS


The New Leica RM2245 – Efficiency Goes Hand in Hand with Comfort

The new Leica RM2245 semi-motorized rotary microtome represents the perfect combination of economic efficiency and user comfort. It is the instrument of choice for clinical histology and histopathology applications, as well as for industrial quality assurance.

The Leica RM2245 is designed to meet the many requirements of modern laboratories. This task is fulfilled without compromise by combining proven technology and exemplary safety standards with functionality and ergonomic design.

Manual sectioning is enhanced by the high-precision motorized specimen feed, which results in efficient operation with maximum section reproducibility. Operator fatigue is reduced through the use of the patented spring system, which eliminates the need for a heavy counter balance weight in the handwheel. Enhanced user comfort results in maximum workplace efficiency and high specimen throughput.

With the Leica RM2245 there is a choice between the conventional manual sectioning method of full handwheel rotation and the “rocking mode”. In rocking mode, the handwheel is only turned back and forth over a short distance. Every change in rotation direction is automatically translated into a specimen feed or retraction movement.



Indiana University Histotechnology Program



The Useful Details Show a Great Design Concept

■ Unobstructed work area

To provide improved access to the cutting area, several clamping levers that are not routinely used have been replaced by hidden clamping devices without compromising standard setting options.

■ Safety during sectioning

Safety is ensured by the integrated, red finger protection guard of the knife holder E that uniquely shields the blade edge. When used correctly the integrated finger protection guard rules out potential injuries because it can even remain on the blade edge during sectioning.

■ Convenient storage area

An integrated storage tray* ensures that ancillary tools are on hand at all times. Alternatively a cooling plate with thermo element can be used for cold storage of specimen blocks either on top or on either side of the microtome.

■ Comfortable operation

All functions are controlled via a compact, space-saving control panel with areas in contrasting colors and touch-sensitive buttons that allow 'blind operation'. Panel inclination is adjustable, and the panel can be positioned either on the left or right side of the instrument.



■ Easy to clean

Due to the rounded shape of the Leica RM2245's one piece housing, section waste is easily removed. The integrated slot cover prevents debris from entering the internal mechanisms of the microtome.



■ Ergonomics and user safety

The smooth-running safety handwheel has an ergonomically shaped handle. A hand-wheel locking device for one-thumb operation locks the specimen head in the uppermost position, which is ideal for changing the specimen and/or blade. A second locking mechanism on the microtome base allows locking the object head in any position.

■ Stability

Maximum rigidity and stability are provided by the optimized microtome base plate.

■ Spacious, magnetized section waste tray

The spacious section waste tray reliably collects sectioning debris. It is safely held in position by magnets, yet is easily removable for emptying and cleaning.

■ Precision specimen orientation with clear zero reference point

This novel orientation mechanism sets new standards of precision and is very easy to operate. An important detail: Calibrated controls make it easy to adjust a specimen to an exact zero position indicated by red pins or a measurable variable on the x/y axis for re-cuts.

■ Lateral knife holder adjustment

The entire width of the microtome blade can be used without having to release the blade clamping mechanism. Using the entire knife edge saves costs, especially in laboratories with high specimen throughput. Three predefined stop positions (left, center, right) correspond to the width of a standard histology cassette.

Well thought-out and proven

- **The Leica-patented Force Compensation System** for extremely smooth manual handwheel operation.
- **The Precision Micrometer Feed System** from 0.5 μm to 100 μm in sectioning mode and 1 μm to 600 μm in trimming mode via stepper motor.
- **The Leica Control Panel Concept** separate and intuitive control panel for instrument settings and communication display integrated in the instrument housing at eye level.
- **Accessories** a wide range of blade or knife holders and specimen holders is available for a variety of applications.

Technical Specifications

Section thickness:

- Setting range: 0.5 – 100 μm
- Setting values: from 0.5 – 5 μm in 0.5 μm -increments
from 5 – 20 μm in 1 μm -increments
from 20 – 60 μm in 5 μm -increments
from 60 – 100 μm in 10 μm -increments
- Trimming section thickness setting range: 1 – 600 μm
- Setting values: from 1 – 10 μm in 1 μm -increments
from 10 – 20 μm in 2 μm -increments
from 20 – 50 μm in 5 μm -increments
from 50 – 100 μm in 10 μm -increments
from 100 – 600 μm in 50 μm -increments

Object feed: 28 mm \pm 1 mm, feed motion via step motor

Vertical specimen stroke: 70 mm

Sectioning modes: 2 manual modes (conventional full handwheel rotation and rocking mode)

Specimen retraction: 5 – 100 μm in 5 μm -increments, can be turned off

Electric coarse feed: 300 $\mu\text{m}/\text{s}$ and 900 $\mu\text{m}/\text{s}$

Maximum specimen size (L x H x W): 50 x 60 x 40 mm

Specimen orientation: horizontal: 8°, vertical: 8°

Nominal supply voltages: 100/120/230/240 V AC \pm 10%

Nominal frequency: 50/60 Hz

Power draw: 70 VA

Dimensions basic instrument:

- Width (including handwheel): 413 mm
- Depth (including waste tray): 618 mm
- Height (with storage area on cover): 305 mm
- Weight (without accessories): approx. 37 kg (87 lbs)

Dimensions control panel:

- W x D x H: 94 x 164 x 50 mm
- Weight: approx. 0.450 kg (1 lb)

Wide range of accessories on request.
Technical specification subject to change.
Certificates: CE, c-CSA-us

Up-to-date development, production and quality control procedures certified under DIN EN ISO 9001 to ensure highest quality and reliability.

The Leica RM2245 features at a glance:

- Precision specimen orientation with zero point reference
- Compact, ergonomic overall design
- Section thickness totalizer and section counter
- Easy alternation between trimming and sectioning function
- Two motorized forward and backward specimen coarse feed speeds
- Programmable specimen retraction
- Low-maintenance cross roller bearings
- Single-handed operation universal cassette clamp
- Completely new knife holder concept for disposable blades, with finger guard in contrasting color
- Precise knife holder lateral movement
- User safety integrated into overall design
- Smooth-running handwheel with integrated quick-lock mechanism
- Ergonomically optimized handwheel handle
- Communication display integrated in instrument housing
- Enclosed micrometer mechanism
- Spacious magnetized section waste tray
- Intuitive control panel
- Wide range of accessories



Leica Microsystems – the brand for outstanding products

Leica Microsystems' mission is to be the world's first-choice provider of innovative solutions to our customers' needs for vision, measurement, lithography and analysis of microstructures.

Leica, the leading brand for microscopes and scientific instruments, developed from five brand names, all with a long tradition: Wild, Leitz, Reichert, Jung and Cambridge Instruments. Yet Leica symbolizes innovation as well as tradition.

Leica Microsystems – an international company with a strong network of customer services

Australia:	Gladesville	Tel. +61 2 9879 9700	Fax +61 2 9817 8358
Austria:	Vienna	Tel. +43 1 486 80 50 0	Fax +43 1 486 80 50 30
Canada:	Richmond Hill/Ontario	Tel. +1 905 762 2000	Fax +1 905 762 8937
Denmark:	Herlev	Tel. +45 4454 0101	Fax +45 4454 0111
France:	Rueil-Malmaison	Tel. +33 1 473 285 85	Fax +33 1 473 285 86
Germany:	Bensheim	Tel. +49 6251 136 0	Fax +49 6251 136 155
Italy:	Milan	Tel. +39 0257 486.1	Fax +39 0257 40 3273
Japan:	Tokyo	Tel. +81 3 5435 9600	Fax +81 3 5435 9615
Korea:	Seoul	Tel. +82 2 514 65 43	Fax +82 2 514 65 48
Netherlands:	Rijswijk	Tel. +31 70 4132 100	Fax +31 70 4132 109
People's Rep. of China:	Hong Kong	Tel. +852 2564 6699	Fax +852 2564 4163
Portugal:	Lisbon	Tel. +351 21 388 9112	Fax +351 21 385 4668
Singapore		Tel. +65 6779 7823	Fax +65 6773 0628
Spain:	Barcelona	Tel. +34 93 494 95 30	Fax +34 93 494 95 32
Sweden:	Sollentuna	Tel. +46 8 625 45 45	Fax +46 8 625 45 10
Switzerland:	Glattbrugg	Tel. +41 1 809 34 34	Fax +41 1 809 34 44
United Kingdom:	Milton Keynes	Tel. +44 1908 246 246	Fax +44 1908 609 992
USA:	Bannockburn/Illinois	Tel. +1 847 405 0123	Fax +1 847 405 0164

and representatives of Leica Microsystems
in more than 100 countries.

The companies of the Leica Microsystems Group operate internationally in four business segments, where we rank with the market leaders.

● Microscopy Systems

Our expertise in microscopy is the basis for all our solutions for visualization, measurement and analysis of microstructures in life sciences and industry. With confocal laser technology and image analysis systems, we provide three-dimensional viewing facilities and offer new solutions for cytogenetics, pathology and materials sciences.

● Specimen Preparation


We provide comprehensive systems and services for clinical histo- and cytopathology applications, biomedical research and industrial quality assurance. Our product range includes instruments, systems and consumables for tissue infiltration and embedding, microtomes and cryostats as well as automated stainers and coverslippers.

● Medical Equipment

Innovative technologies in our surgical microscopes offer new therapeutic approaches in microsurgery.

● Semiconductor Equipment

Our automated, leading-edge measurement and inspection systems and our E-beam lithography systems make us the first choice supplier for semiconductor manufacturers all over the world.

 www.histo-solutions.com

Leica Microsystems Inc.
2345 Waukegan Road
Bannockburn, IL 60015
www.leica-microsystems.com

Telephone 847-405-0123
Toll Free 800-248-0123
Fax 847-405-0164
In Canada call 800-205-3422


MICROSYSTEMS