

Living up to Life



LMD APPLICATION NOTE

Region-specific Gene Expression in Adult Mouse CNS Tissues

Region-specific Gene Expression in Adult Mouse CNS Tissues

Abstract

Different areas of the Central Nervous System (CNS) display specific and selective gene expression profiles. Here, we used the Leica Laser Microdissection system LMD6500 to study region-specific mRNA expression in the adult mouse retina and hippocampus. Quantitative real-time reverse transcription PCR (qRT-PCR) and microarray analyses were performed to verify the specificity of the microdissection procedure and the purity of the extracted RNAs.




Fig.1: LCM of adult mouse retinal layers. A) uncut section, B) microdissection of ONL, C) rhodopsin qRT-PCR from microdissected INL and ONL. Abbreviations: GCL, ganglion cell layer; INL, inner nuclear layer; ONL, outer nuclear layer; Scale bars (A, B): 400 μ m

AUTHORS

Andrea Messina¹, Mark Dunleavy^{2,3}, Paola Sgadò², Valentina Adami^{1,2}, Yuri Bozzi^{2,4},
Simona Casarosa^{1,4}

¹Laboratory of Developmental Neurobiology, Centre for Integrative Biology (CIBIO),
University of Trento, Italy

²Laboratory of Molecular Neuropathology, Centre for Integrative Biology (CIBIO),
University of Trento, Italy

³Department of Physiology, Royal College of Surgeons (RCSI),
Dublin, Ireland

⁴Neuroscience Institute, National Research Council (CNR),
Pisa, Italy

1. Laser capture microdissection of adult mouse retinal layers.

Sections from C57BL/6 adult mouse eyes were cut (15 μ m) with a Leica CM1850UV cryostat and placed on Leica PEN-membrane slides. After hematoxylin staining, LCM was performed with a Leica LMD6500 microdissector using a 5x objective in brightfield with TL-BF contrast method. The laser parameters used for the dissection were: power 60, aperture 7, speed 7, specimen balance 46 and offset 25. Retinal layers were microdissected using the pencil function to border each single layer (Fig. 1A-B). Microdissected layers derived from three mice were pooled, and total RNA was extracted using Macherey-Nagel reagents and protocols. To test the precision of microdissecting single layers, we performed qRT-PCR for rhodopsin, a protein involved in phototransduction specifically expressed of photoreceptors, whose nuclei are located in the outer nuclear layer (ONL) (Fig. 1C).

2. Laser capture microdissection of adult mouse hippocampal subfields.

Sections from C57BL/6 adult mouse dorsal hippocampi were cut (25 μ m) with a Leica CM1850UV cryostat and mounted onto Leica PEN-membrane slides. CA1 and CA3 subregions were dissected using the “Draw and Cut” function of the Leica LMD6500 microdissector (Figure 2A-D). The 10x objective with TL-BF contrast was used. The laser parameters used for the dissection were: power 60, aperture 12,

speed 10, specimen balance 45 and offset 40. To improve the quality and accuracy of cutting, the concave side of the area to be dissected was cut first to prevent the remaining tissue from falling out of the focal plane. The total volume of tissue dissected from each subregion was CA1: 3.1 ± 0.3 mm³, CA3: 3.8 ± 0.4 mm³ (approximately corresponding to 75 sections per brain). Total RNA was extracted from each pool of CA1/CA3 tissue using Qiagen reagents and protocols and used for hybridization of mouse gene expression arrays (Agilent 4X44K). Microarray data analyses by the Agilent GeneSpring GX software and DAVID functional annotational tools showed a significant regulation of genes involved in synaptic function (Fig. 2E-F).

3. Conclusions

Our results confirm that different CNS regions are characterized by specific gene expression profile. LCM, using the Leica system, is an efficient, contact- and contamination-free approach to further investigate CNS differential gene expression, hopefully providing greater insight into a number of CNS pathologies.




Fig.2: LCM of adult mouse hippocampal subfields. A-D) Section of dorsal hippocampus showing the outline and post-cut of the CA3 and CA1 subregions (scale bar: 275 μ m). E-F) Microarray data: bar graph displaying the major biological processes associated with regulated genes (E) and expression heatmap of genes associated with potassium ion transport (F).

The statement by Ernst Leitz in 1907, *“With the User, For the User,”* describes the fruitful collaboration with end users and driving force of innovation at Leica Microsystems. We have developed five brand values to live up to this tradition: Pioneering, High-end Quality, Team Spirit, Dedication to Science, and Continuous Improvement. For us, living up to these values means: **Living up to Life.**

Leica Microsystems operates globally in four divisions, where we rank with the market leaders.

LIFE SCIENCE DIVISION

The Leica Microsystems Life Science Division supports the imaging needs of the scientific community with advanced innovation and technical expertise for the visualization, measurement, and analysis of microstructures. Our strong focus on understanding scientific applications puts Leica Microsystems’ customers at the leading edge of science.

INDUSTRY DIVISION

The Leica Microsystems Industry Division’s focus is to support customers’ pursuit of the highest quality end result. Leica Microsystems provide the best and most innovative imaging systems to see, measure, and analyze the microstructures in routine and research industrial applications, materials science, quality control, forensic science investigation, and educational applications.

BIOSYSTEMS DIVISION

The Leica Microsystems Biosystems Division brings histopathology labs and researchers the highest-quality, most comprehensive product range. From patient to pathologist, the range includes the ideal product for each histology step and high-productivity workflow solutions for the entire lab. With complete histology systems featuring innovative automation and Novocastra™ reagents, Leica Microsystems creates better patient care through rapid turnaround, diagnostic confidence, and close customer collaboration.

MEDICAL DIVISION

The Leica Microsystems Medical Division’s focus is to partner with and support surgeons and their care of patients with the highest-quality, most innovative surgical microscope technology today and into the future.

Leica Microsystems – an international company with a strong network of worldwide customer services:

Active worldwide		Tel.	Fax
Australia · North Ryde	+61	2 8870 3500	2 9878 1055
Austria · Vienna	+43	1 486 80 50 0	1 486 80 50 30
Belgium · Groot Bijgaarden	+32	2 790 98 50	2 790 98 68
Canada · Concord/Ontario	+1	800 248 0123	847 405 0164
Denmark · Ballerup	+45	4454 0101	4454 0111
France · Nanterre Cedex	+33	811 000 664	1 56 05 23 23
Germany · Wetzlar	+49	64 41 29 40 00	64 41 29 41 55
Italy · Milan	+39	02 574 861	02 574 03392
Japan · Tokyo	+81	3 5421 2800	3 5421 2896
Korea · Seoul	+82	2 514 65 43	2 514 65 48
Netherlands · Rijswijk	+31	70 4132 100	70 4132 109
People’s Rep. of China · Hong Kong	+852	2564 6699	2564 4163
· Shanghai	+86	21 6387 6606	21 6387 6698
Portugal · Lisbon	+351	21 388 9112	21 385 4668
Singapore	+65	6779 7823	6773 0628
Spain · Barcelona	+34	93 494 95 30	93 494 95 32
Sweden · Kista	+46	8 625 45 45	8 625 45 10
Switzerland · Heerbrugg	+41	71 726 34 34	71 726 34 44
United Kingdom · Milton Keynes	+44	800 298 2344	1908 246312
USA · Buffalo Grove/Illinois	+1	800 248 0123	847 405 0164