


The World's First Super Zoom Confocal!

Leica TCS LSI



- | | | |
|---|---|--|
| <ul style="list-style-type: none"> Ⓐ Supply unit ① Internal optical bench ② Solid state laser with 405 nm (optional) ③ Laser driver board ④ Power supply ⑤ Acousto-Optical Tunable Filter (AOTF) ⑥ Fiber coupler | <ul style="list-style-type: none"> Ⓑ Scan head ⑦ Confocal excitation pinhole ⑧ Beam splitter, motorized ⑨ Confocal detection pinhole, motorized ⑩ Spectrophotometer prism ⑪ Spectral selector | <ul style="list-style-type: none"> ⑫ Collimation lens ⑬ Photomultiplier ⑭ K-scanner ⑮ Scan lens ⑯ Objective lens ⑰ Transmitted light detector ⑱ Z-zoom lens (motorized, optional) |
|---|---|--|

Leica

MICROSYSTEMS