

システム機器構成セレクションガイド

A - Bの順に選択してください。その他オプション類は必要に応じてお選びください。

A 光学系 (必須)

レーザー光源をVIS(488/532/635nm)または、V-VIS(405/488/532/635nm)から選択してください。

Solid State Laser:

★ R: 635nm 15mW

★ G: 532nm 15mW

★ B: 488nm 15mW

★ V: 405nm 25mW (V-VIS選択時に搭載されます。)

- TCS SPE VIS with HC-Optics
- TCS SPE V-VIS with ACS-Optics

Scanhead

Supply Unit
with integrated PC

RGBV
Solid State AOTF
Laser

Standard Accessories:

Fast Z-Movement
SuperZ Galvo Focus SPE

SuperZ Galvo Inserts
Insert for SuperZ,
standard

Fluo-Illumination:
EL6000

External light source
HX P120 burner
Liquid light guide
Adapter

Input/Output Equipment


TFT Monitor 19"

Keyboard,
Mouse

B 顕微鏡 (必須)


顕微鏡タイプ(正立・倒立)を選択してください。

Microscopes:




DM2500 CSQ VIS
Lasersafe 6x fluo BF
incl. Optical Outfit
and Basic Objective Set

- 正立顕微鏡 DM2500シリーズ
- 15 6601 201 TCS SPE DM2500 VIS
- 15 6601 210 TCS SPE DM2500 V-VIS



DM5500 CSQ VIS
Lasersafe 5x fluo BF
incl. Optical Outfit
and Basic Objective Set

- 正立顕微鏡 DM5500シリーズ
- 15 6601 401 TCS SPE DM5500 VIS
- 15 6601 410 TCS SPE DM5500 V-VIS



DM14000 CSQ VIS
Lasersafe 6x fluo BF
incl. Optical Outfit
and Basic Objective Set

- 倒立顕微鏡 DM14000シリーズ
- 15 6601 501 TCS SPE DM14000 VIS
- 15 6601 510 TCS SPE DM14000 V-VIS

オプション

Additional Options:

Transmitted Light Detectors Kits


- ① 15 6605 100 TLD for DM
- ② 15 6605 101 TLD for DMI

- ③ 15 5931 251 Optical outfit DIC for DM2500 CSQ

- ④ 15 5931 551 Optical outfit DIC for DM5500 CSQ


- ⑤ 15 5933 451 Optical outfit DIC for DM14000 CSQ

Digital Camera Kit



15 5605 300
Leica DFC350 FX incl.
C-Mount Adapter x 0,63(1)

Climate Chamber Kit



15 6605 400
Climate chamber for
DMI4000 CSQ (1)


Including Heating control unit

Including Temp control unit


SuperZ Galvo Inserts

- ⑩ 15 5935 411 Insert for SuperZ, universal
- 15 6604 410 Insert for SuperZ, rotatable
- Option Insert for SuperZ Galvo Focus w/o heating
- Option CO 2 cover for insert for SuperZ Galvo Focus

Tables



15 5934 190
Microscope table for upright and inverted systems, passive.



15 6606 190
Computer table

Note:(1) Available after May 2006

電動ステージ、デジタルカメラの対応時期は別途お問い合わせください。