

The Leica ergonomics program helps people work better

We look after our competitor's customers as well.

We can now offer to users of Nikon and Olympus stereomicroscopes the same comfort that our own customers already enjoy. With the help of the new adapters, you can fit the Leica ErgoTube* 10° - 50° or ErgoWedge* 5° - 25° in place of the existing binocular tube.

*trade mark

Raise your productivity

We have invested more effort than any other stereomicroscope manufacturer to ensure that our instruments are human-engineered. To minimize fatigue resulting from the predominantly sedentary position of the microscopist has long been one of our main goals. With the unequalled

choice of binocular tubes now on offer, the user of a Leica stereomicroscope can match the viewing height and viewing angle to individual requirements. Where the instruments offered by other companies do not offer this facility, the user's position has to be matched to the particular design of the instrument, leading to severe fatigue which often results in lowered productivity.

Invest in ergonomic workstations
Ask us today about our effective ways of improving productivity in science and industry by means of fatigue-free work stations. Enjoy the healthy benefits of the Leica ergonomics program in your daily work.

10 446 251 Tube adapter for Nikon
10 446 250 Tube adapter for Olympus

