

From Eye to Insight

Leica
MICROSYSTEMS

Qualification Services

Installation Qualification (IQ) and
Operational Qualification (OQ) for
your Leica instrument

Be compliant with government
and international standards!



QUALIFICATION SERVICES

Be compliant with government and international standards all over the world!

Who would be better suited to support you with Compliance Services than the experts who develop and support the instruments you use? Our trained and certified service technicians work with you to ensure that your Installation Qualification and Operational Qualification meets the requirements of our overall system validation.

This provides you with well-documented proof that your equipment is installed in accordance with the specifications and meets the specified performance.

Installation Qualification (IQ)

- > The procedure helps ensure that, the instrumentation is installed according to Leica's documented specifications.
- > Completeness of shipping and components of the device is documented.
- > Review and documentation of complete incoming material receipt or compatibility of upgrades with the entire system configuration.



Operational Qualification (OQ)

- > The procedure helps ensure that the instrumentation is operating according to Leica's specifications.
- > Comprehensive testing of the instrument based on defined conditions and known sample characteristics.

Find out more about Qualification Services under
www.leica-microsystems.com/services - and go compliant with us.

Headquarters

Leica Microsystems GmbH
Ernst-Leitz-Strasse 17-37 | D-35578 Wetzlar (Germany) Tel. +49 (0) 6441 29-0 | F +49 (0) 6441 29-2599

www.leica-microsystems.com/service

CONNECT
WITH US!

