

EC Declaration of Conformity
EG Konformitätserklärung



Leica Microsystems CMS GmbH
Ernst-Leitz-Straße 17-37
D-35578 Wetzlar

We hereby declare that the device described below, both in its basic design and construction and in the version marked by us, conforms to the relevant safety- and health- related requirements of the appropriate EU directives.

This declaration shall cease to be valid if modifications are made to the device without our approval.

Hiermit erklären wir, daß nachfolgend bezeichnetes Gerät aufgrund seiner Konzipierung und Bauart sowie in der von uns in Verkehr gebrachten Ausführung, den einschlägigen grundlegenden Sicherheits- und Gesundheitsanforderungen der EU Richtlinien entspricht. Bei einer nicht mit uns abgestimmten Änderung des Gerätes, verliert diese Erklärung ihre Gültigkeit.

Product / Bezeichnung:

Autoloader

Model:/ Gerätetyp:

SL801

Identification No / Gerätenummer :

11544516

EU directives / EU Richtlinien:

In vitro Diagnostic Directive : 98/79/EG

In vitro Diagnostica : 98/79/EG

**Harmonised standards applied /
Angewandte harmonisierte Normen :**

EN 591: 2001
EN 980: 2008
EN 1658: 1996
EN 61000-3-2:2006
EN 61000-3-3:2008
EN 61010-1: 2001
EN 61010-2-081: 2003
EN 61010-2-101: 2003
EN 61326-2-6: 2005
EN 14971: 2007
EN 62304: 2007
ISO 8255-1: 1986
ISO 8037-1: 1986

A handwritten signature in blue ink, appearing to read "Rainer Fuhrmann", with a horizontal line underneath it.

i.V. Rainer Fuhrmann
Qualitätsmanager

Wetzlar, den 25.03.2010