



Product Service

CERTIFICATE

No. Z1 13 01 56257 061

Holder of Certificate: Leica Microsystems CMS GmbHErnst-Leitz-Straße 17-37
35578 Wetzlar
GERMANY**Production Facility(ies):** 60797**Certification Mark:****Product:** Laboratory Equipment
(Microscope)**Model(s):** Leica DMD108**Parameters:**
Microscope:
Rated voltage: 12V
Rated frequency: DC
Rated current: 500mA
Protection class: IIIused with the following external power supplies:
SPU45-105 (SINPRO),
MCEXT12V60W (MULTICOMP) or
MES50A-3 (MeanWell)**Tested according to:** EN 61010-1:2001
EN 61010-2-101:2002
EN 62471:2008

The product was tested on a voluntary basis and complies with the essential requirements. The certification mark shown above can be affixed on the product. It is not permitted to alter the certification mark in any way. In addition the certification holder must not transfer the certificate to third parties. See also notes overleaf.

Test report no.: 028-71320181-300**Date,** 2013-01-09
(Erich Thurner)

Page 1 of 1

